

| | |
|-------------------------------|--|
| Test Name: | Xpert®Xpress CoV-2/ Flu/RSV plus |
| Test Methodology: | Multiplexed real-time RT-PCR |
| Test to be Performed: | Respiratory Multiplex PCR (Flu A/B, SARS-CoV-2, RSV) |
| CPT: | 87637 |
| Reflex Testing Criteria: | If FluA or FluB is positive, sample will be submitted for proper subtyping |
| Pre-Approval Required: | N/A |
| Accepted Specimens: | N/A |
| Specimen Stability: | Specimen should be stored at 2-8°C and tested within 7days of collection. |
| Storage/Transport Conditions: | Specimens should be stored and transported at 2-8°C, or 34-45°F. Sufficient frozen cold packs in an insulated cooler can be used to maintain this temperature. |
| Turnaround Time: | 3 business days |
| Reference Range: | Not detected |
| Comments: | |